

Vein Visualization for Cosmetic Procedures

The AccuVein AV500 is a vein visualization device that displays a map of the vasculature directly on the patient's skin. Locating veins can help with IV starts, blood draw, sclerotherapy, general and cosmetic surgery.



In cosmetic injectable procedures AccuVein helps to:

- Locate and avoid veins to minimize bruising
- Improve patient satisfaction
- Personalize treatment to patient's anatomy

"Seldom does a technological advancement change the way medicine is practiced. In my case, being an Aesthetic Physician injecting hundreds of patients a week, bruising is an unwelcome and dreaded reality. AccuVein has reduced the incidence of bruising to almost zero."

- William H. E. Romero, MD, Cosmetic Physician

In vascular procedures AccuVein helps to:

- Show veins to patients
- Locate feeder veins
- Highlight veins in transdermal laser therapy & sclerotherapy
- Visually confirm efficacy of procedure

"This product changed the way I practice Phlebology. And improved all my sclerotherapy results."

- Dr. Rodrigo Kikuchi, MD, Vascular Surgeon



In PRP procedures AccuVein helps to:

- Identify best access site
- See vein movement or rolling
- Visualize vein patency
- Reduce patient pain and stress

"We have had people who weren't getting a blood draw. We had to send them home. Now with the AccuVein we've been able to draw from them."

- Jina Morys, Technical Director, Acute Care Lab

The Leader in Vein Visualization

AccuVein is the global leader in vein visualization with thousands of facilities around the world using AccuVein's award winning, portable vein visualization. Differentiate your practice and give yourself a competitive advantage with the use of vein visualization by AccuVein.

An Essential Tool

Whether you are performing injections, sclerotherapy, blood draw or any number of surgical procedures, the AccuVein AV500 can help you locate veins so they can either be accessed or avoided. Imagine how much more effective you would be with a map of the vasculature highlighted on the surface of the skin. More effective procedures with less patient discomfort result in higher patient satisfaction, making it clear why so many facilities have chosen to incorporate AccuVein into their standard of care.

AV500 Features

Easy to use	No calibration or adjustment necessary, ready for immediate use. Device can be held in any direction, and works at four or more inches from the skin's surface.
Small size	The AV500 fits in your hand and weighs less than 10 ounces.
Rechargeable battery	No need to remove the battery or plug it into an electrical outlet during use. Charge in cradle or stand and run unplugged for hours.
Viewing modes	In addition to the standard viewing option, the AV500 has an inverse mode which can enhance the visibility of small veins, as well as three brightness settings.
Accuracy	Permanently aligned; center line accurate to less than the width of a human hair.
Patient types	AccuVein's vein illumination will work on a broad range of patients including dark skinned and high BMI patients
Hands-free options	Assess the veins with this lightweight handheld device, then slip it into the "quick cup" of a hands-free accessory freeing up your hands to perform the procedure while viewing the vasculature.

Contact your representative
for product pricing and availability

**Visit www.accuvein.com or call us
for more information at 1.888.909.8888**